Neu ab:
New as of:
Nouveau depuis:
Nuevo a partir de:

07.2012

ORTHOPHOS XG 5

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```
Details.txt x
0 10 20 30 40 50 60 70 80
32
33 -----
34 ----- Standard Configuration DX 11 -----
35 -----
36 Name of flat module: 'DX11 '
37 Code Number: ' 5925214'
38 Serial Number: ' 329039368'
39 Hardware Index: 0
40 Hardware Revision: 2
41 Test Date: 2010-Feb-15
42 Test Counter: 1
43 Processor Derivate: ' MPC855T2P80D4'
44 Previous Software Version: 132
45 Previous Software Revision: 2
46 Date Last Software Update: 2010-Oct-01
47 Current Software Version: 148
48 Current Software Revision: 39
49 GBPS defined area: 0xff 0xff 0xff 0xff 0xff 0xff 0xff 0xff
50
51 -----
52 ----- Extended Configuration DX 11 -----
53 -----
54 DX11 SerialNo: 18757644
55 DX11FpgaVersion: 3
56 DX11FpgaRevision: 4
57 DX11HwVersion: 8
58 DX11MemorySize: 134217728
59 DX11CpuClock: 80020000
60 BspVersion: 1.4/2-IX0204
61
62 -----
63 ----- Error Logging Data DX 11 -----
64 -----
65
66 Timestamp Kategorie Message
67 2010-10-12, 10:13:33 [Message ]: Logbook started
68 2010-10-12, 10:13:55 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
69 2010-10-12, 10:17:22 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:17:22
70 2010-10-12, 10:19:57 [Message ]: Logbook started
71 2010-10-12, 10:20:19 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
72 2010-10-12, 10:21:37 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:18:56
73 2010-10-12, 10:40:27 [Message ]: Logbook started
74 2010-10-12, 10:40:49 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
75 2010-10-12, 10:41:36 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:44:17
76 2010-10-12, 10:59:54 [Message ]: Logbook started
77 2010-10-12, 11:00:15 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
78 2010-10-12, 11:00:41 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 11:00:41
79 2010-10-12, 11:03:10 [Message ]: Logbook started
80 2010-10-12, 11:06:34 [Message ]: Logbook started
```

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Installation Report /
 Warranty Passport

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The undersigned hereby certifies that the in-
 stallation and commissioning of the device
 complies with the requirements of the
 applicable laws and regulations.

Commissioning performed error-free according to
 the installation instructions and installation in-
 structions as well as country-specific safety con-
 ditions where required.
 The unit, including accessories and all docu-
 ments in accordance with the Operating Instruk-
 tions, was turned over to the customer in pri-
 vate and complete condition.
 (USA only: The customer received his copy of the
 Federal Report of Assembly 10-27979.
 The person signing certifies having the required
 technical knowledge for performing commissioning
 and training and fulfilling and complying with the
 requirements outline above.

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07.2012

ORTHOPHOS XG 5

Converting instructions for the "ORTHOPHOS XG 5 to XG 3D upgrade kit"

English


Table of contents

1	About these installation instructions	5
1.1	Scope	5
1.2	Other documentation required.....	5
1.3	Structure of the document	6
1.3.1	Identification of the danger levels.....	6
1.3.2	Formats and symbols used	6
2	Safety instructions	7
2.1	Condensation	7
2.2	Qualifications of service personnel.....	7
2.3	Switching the unit on	7
2.4	Radiation protection	7
2.5	Laser light localizer.....	7
2.6	Electrostatic discharge	8
3	Delivery.....	9
3.1	Scope of supply	9
3.2	3D module	10
3.3	Adjustment sets.....	11
4	Performing a functional test.....	12
4.1	Checking the rotating unit mechanics within a test cycle	12
4.2	Checking the height adjustment mechanics.....	12
4.3	Executing a "Pan symmetry" test exposure	13
5	Procedure for conversion	15
6	Activating the functionality for the Easypad.....	17
7	Mechanical conversion	20
7.1	Preparing the unit for conversion	20
7.2	Removing the cover parts	20
7.3	Removing the headrest and outer ring	23
7.4	Replacing the ring motor	24
7.5	Removing the old diaphragm unit.....	25

7.6	Removing the old X-ray tube assembly	26
7.7	Removing the Pan sensor holder	28
7.8	Attaching the ferrite cores	29
7.9	Installing cable L13 in the plug housing	30
7.10	Laying the light barrier and motor cables on the ring in the sliding plate brackets	31
7.11	Laying fiber-optic cable L15	32
7.12	Attaching cable ties to secure the spiral springs	35
7.13	Installing and connecting the wires for the rotation unit	36
7.14	Installing the 3D module.....	37
7.15	Inserting the 2D sensor	41
7.16	Installing and connecting the wires for the new XG 3D X-ray tube assembly	42
7.16.1	Laying cables on the X-ray tube assembly	43
7.17	Installing and connecting the wires for the new XG 3D diaphragm unit.....	45
7.18	Replacing the user interface	46
7.19	Installing the outer ring and headrest.....	47
7.20	Safety checks.....	47
7.21	Attaching cover parts	47
8	Updating the unit software and configuring the unit	48
9	Adjusting/calibrating the unit	51
10	Final work.....	52
10.1	Attaching the identification and instruction notice	52
10.2	Calling "Extended Details".....	53
10.3	Filling in the installation report and warranty passport.....	55
10.4	Filling in the certificate of conformity	56
10.5	Unit handover	57
11	Disposal/returns	58
11.1	Disposal of removed parts	58
11.2	Returning the X-ray tube assembly to Sirona	58

1 About these installation instructions

1.1 Scope

These installation instructions describe the conversion of the ORTHOPHOS XG 5 / Ceph digital X-ray unit to the ORTHOPHOS XG 3D /Ceph digital volume tomograph. They are intended for use exclusively by trained and authorized distributors and service technicians.

1.2 Other documentation required

In addition to these installation instructions, you need the following documents:

Wiring diagrams

- ORTHOPHOS XG 3D Wiring References: Order No. 63 03 494

Installation Instructions

- GALILEOS / ORTHOPHOS XG 3D Software installation:
Order no. 61 42 389
- GALAXIS Operator's Manual: Order No. 61 23 488
- SIDEXIS XG Digital Radiography Installation Instructions:
Order no. 59 67 356

Service Manual

- ORTHOPHOS XG 3D Service Manual: Order No. 63 03 510

All the documents required are supplied together with the conversion kit on data storage devices (CD/DVD).

1.3 Structure of the document

1.3.1 Identification of the danger levels

To prevent personal injury and material damage, please observe the warning and safety information provided in this document. Such information is highlighted as follows:

DANGER

An imminent danger that could result in serious bodily injury or death.

WARNING

A possibly dangerous situation that could result in serious bodily injury or death.

CAUTION

A possibly dangerous situation that could result in slight bodily injury.

NOTICE

A possibly harmful situation which could lead to damage of the product or an object in its environment.

IMPORTANT

Application instructions and other important information.

Tip: Information on making work easier.

1.3.2 Formats and symbols used

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"> ✓ Prerequisite 1. First action step 2. Second action step or ➤ Alternative action ↪ Result 	Prompts you to do something.
See "Formats and symbols used [→ 6]"	Identifies a reference to another text passage and specifies its page number.
• List	Designates a list.
"Command/menu item"	Identifies commands, menu items or quotations.

2 Safety instructions

2.1 Condensation

Extreme fluctuations of temperature may cause condensation inside the unit. Do not switch the unit on before it has reached normal room temperature. See also Technical Data.

2.2 Qualifications of service personnel

Installation and startup may be carried out only by personnel specifically authorized by Sirona.

2.3 Switching the unit on

Due to the risk of injury caused by malfunction, no person may be positioned in the unit when it is switched on.

2.4 Radiation protection

The valid radiation protection regulations and measures must be observed. The statutory radiation protection equipment must be used.

During an exposure, the service engineer should move as far away from the X-ray tube assembly as the coiled cable of the manual release permits.

With the exception of the service engineer, no other persons are allowed to stay in the room during an exposure.

In case of malfunctions, cancel the exposure immediately by letting go of the exposure release button.

2.5 Laser light localizer

The system incorporates Class 1 laser products.

A minimum distance of 10 cm (4") is required between the eye and the laser. Do not stare into the beam.

Do not use the system with any other lasers, and do not make any changes to settings or processes that are not described in these operating instructions. This may lead to a dangerous exposure to radiation.


2.6 Electrostatic discharge

Electrostatic discharge (abbreviated: ESD – **E**lectro**S**tatic **D**ischarge)

Electrostatic discharge from people can damage electronic components when the components are touched.

Touch a ground point to discharge static electricity before touching any boards.

3 Delivery

3.1 Scope of supply

XG5 upgrade kit (REF 6309467)


- Swim ring
- 3D slide cover, compl.
- Rotation unit, compl.
- 3D module, compl.
- Ring motor, compl.
- Sensor cover
- 3D tube assembly
- 3D diaphragm unit, compl.
- Cable ties (2x)
- Sliding plate
- Bite block foam
- Accessories for fixed patient table
- Hygienic protection
- Geometry phantom, compl. [→ 11]
- XG3D constancy test phantom, compl. [→ 11]
- Ceph adjustment phantom, compl. [→ 11]
- Cables L15, L20, L21 (2x)
- User interface with electronics (Easypad)
- XG3D IT package
- OP XG 3D dongle
- Cable pulling tool (1.5 m rubber hose)
- Technical documents (CD/DVD)
- Tube assembly return note
- Screwdrivers TX10, TX20, TX25

3.2 3D module

Dimensions and weight

Dimensions		Length	Width	Height
	in cm	40	36	23
	in inches	15 3/4	14 1/8	9
Weight	in kg	6		
	in lbs	13 1/4		


Scope of supply


A	3D module
B	NOTICE! The cover plate is only attached loosely.
C	Attachment bracket

3.3 Adjustment sets

3D adjustment/calibration


A	Geometry phantom
---	------------------

C	Ceph adjustment phantom
---	-------------------------

Acceptance and constancy test


A	XG 3D constancy test phantom
---	------------------------------

4 Performing a functional test

4.1 Checking the rotating unit mechanics within a test cycle

The test cycle is executed without radiation. The test cycle is used to check that the unit is functioning correctly and to ensure that a complete, uninterrupted cycle is possible. In the event of increased resistance, the rotating unit stops automatically.


1. Use the arrow keys (A) or (B) to select a PAN program (e.g. P1).
2. Press the T key on the control panel.
↳ The program enters test cycle mode.
3. Press the release button.
↳ The test cycle is started.
4. Check the mechanics of the rotating unit during the test cycle:
↳ Can the cycle be executed without interruption?
↳ Are there any unusual noises during the cycle?
5. If the cycle is interrupted or unusual noises can be heard, the mechanics of the rotating unit are defective.
6. Inform the operator of any deficiencies.
NOTICE! The conversion may only be performed on a unit that is intact.
7. Press the T key on the control panel again.
↳ The program exits test cycle mode.

4.2 Checking the height adjustment mechanics

1. Use the Up/Down keys on the control panel to move the unit up and down through its entire adjustment range.
2. Check the height adjustment mechanics:
↳ Does the height adjustment function correctly?
↳ Are there any unusual noises during the adjustment?
3. If the height adjustment does not function correctly or unusual noises can be heard during the adjustment (e.g. a hammering noise which changes with the speed), the mechanics of the height adjustment or the height adjustment motor are defective.
4. Inform the operator of any deficiencies.
NOTICE! The conversion may only be performed on a unit that is intact.

4.3 Executing a "Pan symmetry" test exposure

1. Execute a "Pan symmetry" test exposure (refer to the ORTHOPHOS XG 3D Service Manual, "Maintenance" → "Checking the diaphragm" → "Pan symmetry test exposure").

A	Adjustment OK	<ul style="list-style-type: none"> • The shadow of the center needle, the needle image, and the auxiliary lines must be coincident and located one behind the other. A tolerance (offset of needle from the central auxiliary line) of ± 0.75 mm is admissible. • Distance A1 must be 88.6 ± 1 mm. • Distances A2 must be identical, each being 44.3 ± 0.5 mm.
B		<ul style="list-style-type: none"> • A surrounding white border must be present.
C	Adjustment not OK	

2. Examine the exposure.
3. If there is something wrong with the exposure, inform the operator.
NOTICE! The conversion may only be performed on a unit that is intact.

5 Procedure for conversion

IMPORTANT

The unit may only be converted once a function check [→ 12] has been performed and no deficiencies have been found on the unit.

Tip: Sirona recommends calling [→ 53] the "Extended details" before starting the conversion in order to read out the current parameters from the unit and storing these as an XML file under the network name of the unit in the PDATA/.../P2K_Config directory for documentation purposes.

Activation key code

- Have the activation key code ready for the conversion.
The activation key code can be requested from the Sirona Customer Service Center (CSC).

Activating the functionality for the Easypad

- Activate the functionality for later use of the Easypad [→ 17].

Installing the IT package

- Install the IT package (see GALILEOS / ORTHOPHOS XG 3D Software Installation technical documentation (REF 61 42 389)). This can be found on the documentation CD in the IT package.

Mechanical conversion

1. Remove the Pan sensor and move the unit to a comfortable working position [→ 20].
2. Remove the cover parts [→ 20].
3. Remove the head support and outer ring [→ 23].
4. Replace the ring motor [→ 24].
5. Remove the old diaphragm unit [→ 25].
6. Remove the old X-ray tube assembly [→ 26].
7. Remove the Pan sensor holder [→ 28].
8. Install the ferrite cores [→ 29].
9. Install the connection socket for cable **L13** in the housing [→ 30].
10. *If not yet laid:*
Lay light barrier cable **L21** and motor cable **L20** on the ring in the brackets of the sliding plate [→ 31].
11. *If not yet laid:*
Lay [→ 32] fiber-optic cable **L15**.
12. Attach cable ties to secure the spiral springs [→ 35].
13. Install and connect the wires for the rotation unit [→ 36].
14. Install and connect the wires for the 3D module.
15. Insert the 2D sensor.
16. Install and connect the wires for the new XG 3D X-ray tube assembly [→ 42].
17. Install and connect the wires for the diaphragm unit.
18. Replace the user interface [→ 46].
19. Install the outer ring and head support [→ 47].
20. Attach the cover parts [→ 47].

Updating the unit software and configuring the unit

- Update the unit software and configure the unit [→ 48].

Adjusting the light localizer

1. Adjust the FH (Pan) double laser light localizer (see ORTHOPHOS XG 3D Service Manual, "Repairs → Laser light localizers" chapter).
2. Adjust the MS (Pan) laser localizer (see ORTHOPHOS XG 3D Service Manual, "Repairs → Laser light localizers" chapter).

Adjusting/calibrating the unit

- Perform a complete unit adjustment or calibration [→ 51].

Final work

- Complete the unit conversion as described in the chapter entitled "Final work [→ 52]" .

Disposal of removed parts




- Dispose of parts which are no longer required as described in the chapter entitled "Disposal/returns [→ 58]" .

6 Activating the functionality for the Easypad

Perform a software update

1. Switch the unit on.
2. In *"Automatic"* update mode, perform an automatic software update to the latest XG3D unit software version supplied, as described in the ORTHOPHOS XG5/XG^{Plus} and XG3D^{ready} Service Manual.
NOTICE! At this point, do not acknowledge any error messages on the unit or the PC (e.g. "Timeout suspected").
3. Switch the unit off.
Wait approx. 1 minute. Then switch the unit back on.
Acknowledge any error messages that appear with the R key.


Activating the service menu


1. Press and hold down Service key (D) until the patient symbol keys light up (E-H) (approx. 2 s).
2. Then press the patient symbol keys in the sequence F – H – E within the next 4 s.

➡ After you have entered the key combination correctly, the service menu is displayed.

Activating and executing service routine S008.30 (with security access)


1. Select service routine S008 in selection field 1 (A) with the arrow keys and confirm your selection by pressing the Service key (B).



2. Select test step 30 in selection field 2 with the arrow keys (C) and confirm your selection by pressing Service key (B).

➔ After you have confirmed your selection, a "0" appears in selection field 2.

3. Confirm security access by once again selecting the number of the main routine ("8") in selection field 2 (C) with the arrow keys and press the Service key (B) to confirm your selection.

↗ The message "Irreversible!" flashes on the display.

↗ The LED above the Memory key (D) lights up.

4. Press the Memory key (D).

↗ The LED above the R key lights up.

5. Press the R key (E).

↗ The message "OK" appears on the display.

6. Switch off the unit.

7. Now perform the mechanical conversion of the unit [→ 20].

7 Mechanical conversion

7.1 Preparing the unit for conversion

1. Switch the unit on and move it to a comfortable working position (board DX1 must be easily accessible).
2. Switch the unit off again.

WARNING


Potentially lethal shock hazard!

It is essential to switch the unit off and to wait at least another 4 minutes before removing a cover.

3. Remove the Pan sensor.
Keep this in a safe place for reassembly.

7.2 Removing the cover parts

When removing covers, always remember that direct sunlight or bright room lighting can cause system malfunctions due to activated light barriers. Therefore: avoid direct sunlight and bright room lighting above the unit!


- Remove the cover parts as shown in the drawing.
NOTICE! Before removing the "slide cover, complete", detach connector X607 of the light localizer cable from board DX1.


A	Arm cover, top
B	Slide cover, top rear
C	Slide cover, bottom rear
D	Slide cover, complete (is replaced by the 3D slide cover (included) during reassembly.)
E	Sensor holder cover, large
F	Ring cover
G	Tube assembly cover, front
H	Tube assembly cover, rear
I	Support cover, bottom

IMPORTANT: The cover parts "slide, complete" and "sensor holder, large" are no longer required for reassembly. These are replaced by new cover parts included in the conversion kit.
Dispose of parts which are no longer required in accordance with the regulations applicable in your country [→ 58].

7.3 Removing the headrest and outer ring


1. **NOTICE! Cable L18**
While holding the headrest firmly from below, loosen the screw (A) and remove the headrest.
2. Pull cable **L18** off connector **X1** on board **DX5**.

3. Loosen the four screws (B) and remove the outer ring.
4. Keep the headrest and the outer ring in a safe place for reassembly.



7.4 Replacing the ring motor

1. Detach the ring motor cable from the cable harness and pull it off connector **X813** on board **DX1**.
2. Loosen the four screws (**A**) on the ring motor and remove the motor including the screws and the serrated washers (**B**).
3. Insert the new motor into the ring.
Tip: While inserting the motor, turn it back and forth slightly until the pinion engages in the ring gear.
4. Use the screws (**A**) and serrated washers (**B**) to screw the new motor onto the ring motor support securely.
5. Run the ring motor cable along its original path and plug it back into connector **X813** on board **DX1** (see the ORTHOPHOS XG 3D/Ceph Service Manual).
IMPORTANT: Don't forget to reattach all cable ties and clamps.

7.5 Removing the old diaphragm unit

1. Pull cables **L11** and **L12** off connectors **X501** (L11) and **X101** (L12) on board **DX61**.
2. Loosen the upper screw (**A**) and the two lower screws (**B**) and remove the diaphragm unit.


IMPORTANT

If the two lower screws are not accessible, push the lower slider (**C**) upward manually. If the slider is in such an unfavorable position that it cannot be pushed upward:

- remove board **DX61**,
- loosen the four screws on the motor (**D**) and raise it with the spindle.

NOTE: The removed diaphragm unit is no longer required for reassembly. Dispose of the old diaphragm unit in accordance with the regulations applicable in your country [→ 58].


7.6 Removing the old X-ray tube assembly

1. Loosen the four screws (A) and remove cover plate (B) incl. the cable shielding (L3). **CAUTION! Also pull cable L3 off connector X3 and the ground cable off connector X304 on board DX6.**

Tip: The ferrite core and cable shielding can remain on the cover plate.

2. Detach cables L5 and L6 from the rubber grommets and pull the cables out of sockets J6 (L5) and J2-J3 (L6) on board DX6.

3. Loosen the two rear screws (C) on the tube assembly.

4. CAUTION! The tube assembly is heavy!


Hold the tube assembly firmly in place, loosen the two front screws (D) (3-4 turns) and remove the tube assembly toward the front.

IMPORTANT


The old X-ray tube assembly must not be reused.

Within the scope of the upgrade, it must be packed into the packaging of the new X-ray tube assembly and returned to Sirona with the enclosed return note (see chapter entitled Returning the X-ray tube assembly to Sirona [→ 58]).

7.7 Removing the Pan sensor holder

1. **CAUTION! Cable L13!**
Loosen the three screws (C) and remove the sensor holder.
2. Loosen the two screws (D) and remove the "small sensor holder" cover.
3. Loosen the three screws (E) and remove the connection socket.
4. Remove cable L13 from the shield terminal (F) at the rear of the connection socket.

5. Loosen the connector screw connection on the front side of the connection socket and thread the connector to the rear through the connection socket and the sensor holder.

NOTE: The removed sensor holder is no longer required for reassembly. Dispose of the sensor holder in accordance with the regulations applicable in your country [→ 58].

7.8 Attaching the ferrite cores

DANGER

Potentially lethal shock hazard!


Switch off the power supply connection before the power cable (D) is disconnected or connected.


DANGER

Potentially lethal shock hazard!

Be sure to connect the second protective ground wire (C).


1. Fit the HA motor ferrite core (A) to the motor cable and attach this with a cable tie to the sheet metal as shown in the figure to the left.
2. Loosen the three cores of the power cable (L = brown, N = blue, PE = green/yellow) from the right side of the power line terminal. From underneath, thread the cores twice through the ferrite core mains filter (B) (in a loop).
3. Loosen the second, additional protective ground wire (C) (green/yellow) from the stand.
4. Feed the second protective ground wire through the ferrite core (B) from underneath.
5. Secure the cores with a cable tie to the ferrite core (B) as shown in the drawing (right).
6. Connect the cores of the power cable (L, N and PE) to the power line terminal once more.
7. Re-attach the second protective ground wire (C).


C	Second protective ground wire
D	Power cable to the unit
E	Unit

7.9 Installing cable L13 in the plug housing

1. **NOTICE!** Ensure that the socket is aligned correctly in the housing. If the cable outlet points to the bottom right (see photo), the wider side of the socket (A) must point forwards.
Insert cable L13 into the underside of the connection socket.
2. **NOTICE!** Ensure that both phase elements of cable L13 are secured by the strain relief.
Attach the terminal (B) as the strain relief in the housing.
3. Close the plug housing.
4. Screw the hexagon head bolts (for securing the connector) into the plug housing.

7.10 Laying the light barrier and motor cables on the ring in the sliding plate brackets


Condition before conversion

Replacing old cable fixings with sliding plate

1. Remove **all** old (i.e. already present on the unit) cable fixings (A).

2. Insert the sliding plate (B) behind the edge of the rotary ring, push it up, and clip it (C) to the rotary ring at both sides.


Condition following conversion

Installing cables L20 and L21


3. Use the two screws (D) to screw the sliding plate firmly onto the rotary ring.
 4. Remove the cables from the old cable fixings and lay them in the brackets (E) of the sliding plate in the same order as before.
- *Only if not yet laid:*
Lay cables L20 (1x) and L21 (2x) (included in the upgrade kit) in the sliding plate brackets.

7.11 Laying fiber-optic cable L15

Pulling fiber-optic cable L15 through the cable conduit of the rotating unit


1. Release the cable conduit which accommodates cable **L13** on one side from the rotating unit. This makes fiber-optic cable **L15** easier to install (see figure below).
2. Secure fiber-optic cable **L15** to the cable pulling tool (rubber hose included in the upgrade kit).
Tip: Make a small cut at the end of the rubber hose, push the connector for fiber-optic cable **L15** into the hose and secure the connection with adhesive tape.


3. Carefully feed fiber-optic cable **L15** into the cable conduit (which also accommodates cable L13) from above using the cable pulling tool.
4. Release fiber-optic cable **L15** from the cable pulling tool.
5. Re-attach the cable conduit of the rotating unit.

Laying and connecting fiber-optic cable L15 (DX1 side) from the rotating unit in the arm to board DX1


1. Guide fiber-optic cable **L15** (DX1 side) from below through the cable conduit upward into the arm and route it parallel to the other fiber-optic cables under the cable attachments to board **DX1**.
2. Attach the connector for fiber-optic cable **L15** to socket **J309** on board **DX1**.


Laying fiber-optic cable L15 (DX6 side) on the ring

3. Attach a radius limiter (included in the upgrade kit) level with the cable feed-through in the arm (A) and in front of the plug connection to board DX1 (B) on fiber-optic cable L15.
- Route fiber-optic cable L15 (DX6 side) parallel to cables L20 and L21 along the ring in the sliding plate brackets to the tube assembly.
- IMPORTANT:** The cables should be laid in the brackets in the following order: from below: L20, L21, L15, L21. The thinner fiber-optic cable L15 should be second from the top so that it cannot slip out of the sliding plate brackets.

7.12 Attaching cable ties to secure the spiral springs

1. Check whether the spiral springs in the two corrugated tubes are secured with cable ties (A) at the sensor and tube assembly sides to prevent them slipping out.

NOTICE

The closures for the two cable ties must be on the outer side of the relevant corrugated tube so that they do not damage the cables on the inside if the tube moves.

2. If this is not the case, use an awl to pierce the corrugated tubes at the indicated position (B) (working from the inside out) and fasten each spring to the relevant corrugated tube using a cable tie.



7.13 Installing and connecting the wires for the rotation unit

CAUTION

The unit is rotated via a motor drive.

Never turn the rotation unit manually, as this could damage the gearing.

1. Screw the rotation unit onto the ring using the three screws (A) and washers (B) and the connector plate (C).

2. Plug cable L20 and both light barrier cables L21 (from the rotation unit) into connectors X203 (L20), X303 (L21) and X304 (L21) on the ring.

3. Plug cables L13 and L13a together through the connector plate (C) and screw the plug connection tight.

7.14 Installing the 3D module

NOTICE


Damage to the flat panel

The flat panel of the 3D module is a sensitive component. Be careful not to touch the flat panel during assembly.


NOTICE

Damage to the rotation unit


The 3D module is rotated via a motor drive. The gearing of the 3D module can be damaged if it is turned by hand.


1. Carefully remove the elements of the 3D module from the packaging.
2. **NOTICE!** The EMC cover is only attached to the 3D module on delivery. Remove the EMC cover (A) from the 3D module.


3. Attach the "top 3D module" cover:
 Loosen the screw (B) (approx. 2-3 turns).
 Thread the ring of the "top 3D module" cover over the rotation unit and slide the cover under the ring cover.

4. Screw the cover tight to the rotation unit with the two screws (**K**) (included with the installation material).
Re-tighten screw **B**.

5. Hook the 3D module (**C**) onto the rotation unit (**D**) and secure it immediately with the 3 screws (**E**).


6. Attach the sheet metal bracket (**F**) (incl. cables and strain reliefs) to the 3D module using the 3 screws (**G**).


X212
X211
X210
X209


7. Plug the cables into the sockets of board **DX88** according to their color coding:
X212 = white
X211 = black
X210 = red
X209 = blue

8. Reattach the EMC cover (A) and fasten it with the 4 screws (H).

9. Screw the retaining bracket (I) onto the 3D module with the 4 screws (J).


7.15 Inserting the 2D sensor

NOTICE


Damage to the rotation unit

Insert the 2D sensor carefully into the sensor holder from below and avoid any unnecessary leverage.


- **NOTICE!** The lock must engage audibly.
Insert the removed 2D (Pan) sensor (K) into the sensor holder of the 3D module.

7.16 Installing and connecting the wires for the new XG 3D X-ray tube assembly


1. Hang the new tube assembly on the two front screws (D) of the rotating unit and tighten them securely.
2. Insert the two rear screws (C) and tighten them firmly.
3. Plug cables L3, L5 and L6 as well as the ground cable back onto board DX6 and reattach the cables to the rubber grommets.


4. Re-attach the cover plate (B) with the four screws (A).

7.16.1 Laying cables on the X-ray tube assembly


➤ Lay cables L20, L21 and L15 as shown in the photos below.


A	Cable L5 → Socket J6 on board DX6
B	Cable L6 → Socket J2/J3 on board DX6
C	Grounding cable → Connector X304 on board DX6
D	Cable L3 → Connector X3 on board DX6
E	Laying cables correctly on the cover plate
F	Cable L15 → Socket J5 on board DX6
G	Cable routed on left side of tube assembly: 2x L21, L20 and L11.
H	Cable L12 routed on right side of tube assembly.

7.17 Installing and connecting the wires for the new XG 3D diaphragm unit


1. **NOTICE!** Ensure that the new diaphragm unit is securely hooked into the tube assembly with the bottom two flaps.
Hook the new diaphragm unit with the bottom two flaps into the brass rings (B) on the tube assembly and fix them again with the screw (A).
2. Plug cables L20 and L21 into connectors X203 (L20), X302 (L21) and X303 (L21) on board DX61.
3. Plug cables L11 and L12 into connectors X501 (L11) and X101 (L12) on board DX6.


7.18 Replacing the user interface

B	Easypad
C	Multipad

1. Press into the slit (A) of the cover with a screwdriver (do not pry!) and remove the Multipad from the control panel.
2. Pull cables L9 and L10 off connectors X102 (L9) and X103 (L10) on board DX71 of the Multipad.
3. Plug the cables into connectors X102 (L9) and X103 (L10) of board DX7 on the Easypad.
4. Clip the Easypad onto the control panel.
5. Update the nameplate at the control panel cover.
To do so, affix the supplied label as shown in the figure.

6. *If this does not already exist:*
 Attach a ferrite core (D) (included in the scope of supply for the control panel) to cable L10 (green cable).

7.19 Installing the outer ring and headrest

1. Attach the outer ring incl. swim ring [→ 23].
2. Assemble the headrest [→ 23].

7.20 Safety checks

Once the installation/conversion work is complete, the protective ground conductor and unit leakage current must be tested (see Service Manual, "Maintenance" chapter).

7.21 Attaching cover parts

IMPORTANT! Instead of the old slide cover, the "slide cover, compl.", install the new "3D slide cover, compl." (included in the scope of delivery).

- Install the cover parts [→ 20].

8 Updating the unit software and configuring the unit

IMPORTANT: The current update file can be found on the ORTHOPHOS XG3D Firmware CD. This is included in the XG 3D conversion set (IT package) and also the country-specific set. The contents of the CD can be downloaded from the Dealer domain of the Sirona Internet home page (under Product Info/X-ray Systems): www.sirona.com

Perform a software update


- ✓ The functionality for use of the Easypad must be activated [→ 17].
 - ✓ The IT package from the XG3D upgrade kit must be installed.
 - ✓ The unit is mechanically converted to ORTHOPHOS XG 3D.
1. In *"Automatic"* update mode, perform an automatic software update to the latest unit software version, as described in the ORTHOPHOS XG5/XG^{Plus} and XG 3D^{ready} Service Manual.

NOTICE! At this point, do not acknowledge any error messages on the unit or the PC (e.g. "Timeout suspected").


2. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.

✎ Error message **E6 11 07** (undefined system class) is displayed.

Confirming the unit class

1. Acknowledge the error message with the R key .
✎ The access level for the service menu (level 4) is automatically started.
2. Press and hold down the service key (D) until the patient symbol keys light up (approx. 2 s).
3. Then press the patient symbol keys in the sequence **F – H – E** within the next 4 s.

Confirming the unit serial number

Once the key combination has been entered correctly, service routine S017.1 (select/confirm unit class) is started automatically.

The memory key (A) lights up.

4. Confirm the unit class "01":
To do so, touch the memory key (A) (R key on the touchscreen lights up) and then the R key on the Touchscreen (B).
 5. Exit the service routine with the double-arrow key (C).
 6. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.
- Error message **E61505** (undefined unit serial number) is displayed.

1. Acknowledge the error message using the R key on the control panel.
2. Call the Service menu
(see ORTHOPHOS XG 3D Service Manual).
3. Call service routine S008.3 (see ORTHOPHOS XG 3D Service Manual, "Service → S008: Test step 3").
4. Check the serial number
IMPORTANT: The unit serial number is located on the nameplate of the unit.
NOTICE! If the serial number is incorrect, cancel the update and contact the Sirona Customer Service Center!
5. Confirm the serial number using service routine S008.3 (see ORTHOPHOS XG 3D Service Manual, "Service → Basic operating sequences in the Service menu").
6. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.

Defining the unit version

1. Call service routine S017 and set the unit version using service routines S017.2 (see ORTHOPHOS XG 3D Service Manual, "Service → S017: Test step 2").
 2. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.
- The message "No Key" may appear on the Easypad.

Entering the activation key code

1. Enter the activation key code using the Easypad (see ORTHOPHOS XG 3D Operating Instructions, "Settings and repairs → Activate function → Activation via Easypad").
 2. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.
- The message "No Key" must no longer appear.

Selecting a language and language set

- Call service routine S017 and set the language and the language set using service routines S017.4 and S017.5 (see ORTHOPHOS XG 3D Service Manual, "Service → S017: Test step 4 and 5").

Adding a unit as an XG 3D X-ray component

1. Call the SİDEXİS Manager, delete the XG X-ray component and add the unit as an XG 3D X-ray component.
2. **NOTICE! Now the unit is identified as ORTHOPHOS XG 3D.**
Repeat the software update to the current unit software version (in update mode "*Automatic*"). Proceed as described in the ORTHOPHOS XG 3D Service Manual, "General operating sequences → Software update".

Saving data from board DX61 to board DX11

- Call service routine S017.17 and transfer the data from board DX61 to board DX11 (see ORTHOPHOS XG 3D Service Manual, "Service → S017: Test step 17", Case 1).
IMPORTANT: The data direction shown in selection field 3 must be **DX61 → DX11**.

Saving data to board DX88


1. Call service routine S009.8 (saving/restoring DX88 data) and execute the routine (see ORTHOPHOS XG 3D Service Manual, "Service → S009: Test step 8").
 2. Switch the unit off.
Wait approx. 1 minute. Then switch the unit on again.
- 👉 The unit software has now been updated and the unit is configured. Following a unit restart, only error message E1 11 20 (unit is not calibrated) should appear.

9 Adjusting/calibrating the unit

1. In SIDEXIS XG, call the "Adjustment/calibration" menu:
*"Utilities" / "Constancy test..." / "3D" / ("Select X-ray device") /
 "Service exposure" / ("Select X-ray component") / "Adjustment/
 calibration" / Password request (see also Service Manual, chapter
 "Adjusting/calibrating the unit")*

IMPORTANT

Before you begin adjusting or calibrating the unit, you must reset the 3D adjustments.


2. In the structure tree, under "3D Reset adjustment", click on the "Reset adjustment" element.
3. Click on the "Reset values" button.
 - ↳ **IMPORTANT!** After resetting the values, the parameter **C1** must have a value of "88".
4. Now perform a complete unit adjustment/calibration, as described in the ORTHOPHOS XG 3D Service Manual, "Unit adjustment/calibration" chapter.

10 Final work


10.1 Attaching the identification and instruction notice

- Attach information labels to the covers of the unit. This involves gluing the labels supplied in position as shown in the illustrations.

On the cover of the X-ray tube assembly


On the stand

10.2 Calling "Extended Details"

1. Start the *"SIDEXIS Manager"* under "Start"/"Programs"/"SIDEXIS"/*"SIDEXIS XG"*.
2. Click on *"Configuration of the X-ray components"*.
 - ↳ The *"Configuration of the X-ray components"* menu opens.

3. Select the *"Attributes"* tab.
4. Click on the *"Details..."* button.
 - ↳ The current parameters are read from the unit and filed in an XML file under the network name of the unit in the PDATA/.../P2K_Config folder. The process can take up to 30 seconds. After the parameters are read, an editor displaying the data is opened automatically.

```

Details.txt x
0 10 20 30 40 50 60 70 80
32 -----
33 ----- Standard Configuration DX 11 -----
34 -----
35 -----
36 Name of flat module: 'DX11 '
37 Code Number: ' 5925214'
38 Serial Number: ' 329039368'
39 Hardware Index: 0
40 Hardware Revision: 2
41 Test Date: 2010-Feb-15
42 Test Counter: 1
43 Processor Derivate: ' MPC855T2P80D4'
44 Previous Software Version: 132
45 Previous Software Revision: 2
46 Date Last Software Update: 2010-Oct-01
47 Current Software Version: 148
48 Current Software Revision: 39
49 GBPS defined area: 0xff 0xff 0xff 0xff 0xff 0xff 0xff 0xff
50
51 -----
52 ----- Extended Configuration DX 11 -----
53 -----
54 DX11 SerialNo: 18757644
55 DX11FpgaVersion: 3
56 DX11FpgaRevision: 4
57 DX11HwVersion: 8
58 DX11MemorySize: 134217728
59 DX11CpuClock: 80020000
60 BspVersion: 1.4/2-IX0204
61
62 -----
63 ----- Error Logging Data DX 11 -----
64 -----
65
66 Timestamp Category Message
67 2010-10-12, 10:13:33 [Message ]: Logbook started
68 2010-10-12, 10:13:55 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
69 2010-10-12, 10:17:22 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:17:22
70 2010-10-12, 10:19:57 [Message ]: Logbook started
71 2010-10-12, 10:20:19 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
72 2010-10-12, 10:21:37 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:18:56
73 2010-10-12, 10:40:27 [Message ]: Logbook started
74 2010-10-12, 10:40:49 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
75 2010-10-12, 10:41:36 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:44:17
76 2010-10-12, 10:59:54 [Message ]: Logbook started
77 2010-10-12, 11:00:15 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
78 2010-10-12, 11:00:41 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 11:00:41
79 2010-10-12, 11:03:10 [Message ]: Logbook started
80 2010-10-12, 11:06:34 [Message ]: Logbook started

```

10.3 Filling in the installation report and warranty passport

IMPORTANT

Please note that any warranty claims can be fulfilled only if the fully and legibly completed "Installation Report / Warranty Passport" document is available to Sirona immediately after installation.

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1. Enter the serial numbers and software version in the attached document "Installation Report and Warranty Passport".
2. Complete the remaining points together with the customer.
3. Hand over the copy marked "For the customer" to the customer.
IMPORTANT: In Germany, the copy marked "For the customer" must be kept in the X-ray system book.
4. Fax the "For the dealer" copy back to Sirona immediately after the installation.
IMPORTANT: The copy marked "For the dealer" is for the customer records kept by the dealer.

10.4 Filling in the certificate of conformity

**ORTHOPHOS XG 3D / Ceph,
ORTHOPHOS XG 3D^{ready} / Ceph**

Compatibility list / Declaration of conformity by system integrator
for digital radiography with SIDEXIS

sirona
The Dental Company

System integrator (system installer)		Operator	
Company, address		Name, address	
Component			
PC unit	Manual	Requirement / Approval	Technical characteristics
	Type	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	IBM compatible PC with DualCore 2 GHz RAM ≥ 2 GB Windows XP Professional / Ultimate ≥ 500GB free storage space Graphics card ≥ 1280 x 1024, ≥ 16,7 Mbit colors (True Color), ≥ 512 MB Refresh rate > 70Hz Operating system: according to the op. systems approved in the installation instructions. Removable disk drive¹: DVD-RW Dual Layer
	No		
Monitor	Manual	Requirement / Approval	Technical characteristics
	Type	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	Min. screen diagonal 15" with full screen, 17" with CRT Min. resolution 1280 x 1024 Refresh rate > 70Hz CRT: Dot mask min. 0,28mm flat screen monitor: Pixel pitch 0,30 mm x 0,30 mm
	No		
ORTHOPHOS XG 3D / Ceph, ORTHOPHOS XG 3D^{ready} / Ceph	Manual	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC
	Type		see installation instructions
	No		
	Manual	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC
	Type		see installation instructions
	No		
Other devices optionally connected to the PC (printer, LAN...)	Manual	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC
	Type		see installation instructions
	No		
	Manual	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC
	Type		see installation instructions
	No		
Other devices optionally connected to the PC (printer, LAN...)	Manual	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC
	Type		
	No		
	Manual	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC
Type			
No			
Manual	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC
Type			
No			

1. supported by SIDEXIS, function required only if not performed by an additional server

Declaration acc. to Medical Devices Directive (MDD) 93/42/EEC and IEC 60 601-1-1 item 3.201.4

The undersigned herewith certifies that the above system consisting of CE marked medical devices has been assembled in accordance with the present manufacturer's instructions, that it has been tested for its intended use and that all relevant user information including pertinent information provided by the corresponding manufacturers has been provided for placing it on the market, and that this declaration shall be kept on hand for the competent authorities for a period of five years.

Place, date

Name of system integrator (in block letters)

Signature

D 3352.128.03.01.02

10.2010

A. Nr.: 000 000

63 03 593 D 3352

Operator's copy

- Fill in the certificate of conformity.
The system integrator uses the compatibility list as evidence of the installed SIDEXIS XG system's EC-conformity.
The list is structured in such a way that only the white fields need to be completed. The original list remains in the practice, the copy is given to the system integrator.
- Perform the acceptance test prescribed in the X-ray Ordinance §16. Use the test phantom, plug-in plate, contrast elements, needle phantom and relevant instructions included in delivery for this purpose.

Only for Germany

10.5 Unit handover

- Give the customer the technical documentation as well as all of the patient positioning aids, test phantoms and special tools included in delivery, including packaging. These are important components of the unit and must be stored carefully.

11 Disposal/returns

11.1 Disposal of removed parts

Disposal of removed parts

- Ensure that the removed parts are subjected to environmentally sound recovery and disposal. These parts must not be disposed of as domestic refuse! Please observe the disposal regulations applicable in your country.

11.2 Returning the X-ray tube assembly to Sirona

1. Complete the return note included in the conversion kit.
2. **NOTICE! Use the packaging material in which the new X-ray tube assembly from the conversion kit was delivered.**
Package the old X-ray tube assembly.
Include the completed return note in the packaging.
3. Return the X-ray tube assembly to the Sirona Returns Service in Bensheim.

We reserve the right to make any alterations which may be required due to technical improvements.

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```

Details.txt x
0 10 20 30 40 50 60 70 80
32 -----
33 ----- Standard Configuration DX 11 -----
34 -----
35 -----
36 Name of flat module: 'DX11 '
37 Code Number: ' 5925214'
38 Serial Number: ' 329039368'
39 Hardware Index: 0
40 Hardware Revision: 2
41 Test Date: 2010-Feb-15
42 Test Counter: 1
43 Processor Derivate: ' MPC855T2P80D4'
44 Previous Software Version: 132
45 Previous Software Revision: 2
46 Date Last Software Update: 2010-Oct-01
47 Current Software Version: 148
48 Current Software Revision: 39
49 GBPS defined area: 0xff 0xff 0xff 0xff 0xff 0xff 0xff 0xff
50
51 -----
52 ----- Extended Configuration DX 11 -----
53 -----
54 DX11 SerialNo: 18757644
55 DX11FpgaVersion: 3
56 DX11FpgaRevision: 4
57 DX11HwVersion: 8
58 DX11MemorySize: 134217728
59 DX11CpuClock: 80020000
60 BspVersion: 1.4/2-IX0204
61
62 -----
63 ----- Error Logging Data DX 11 -----
64 -----
65
66 Timestamp Categorie Message
67 2010-10-12, 10:13:33 [Message ]: Logbook started
68 2010-10-12, 10:13:55 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
69 2010-10-12, 10:17:22 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:17:22
70 2010-10-12, 10:19:57 [Message ]: Logbook started
71 2010-10-12, 10:20:19 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
72 2010-10-12, 10:21:37 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:18:56
73 2010-10-12, 10:40:27 [Message ]: Logbook started
74 2010-10-12, 10:40:49 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
75 2010-10-12, 10:41:36 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 10:44:17
76 2010-10-12, 10:59:54 [Message ]: Logbook started
77 2010-10-12, 11:00:15 [Error ]: E1 11 20 (ERR_DX11_INVALID_CB_CALIBDATA)
78 2010-10-12, 11:00:41 [RTC Date/Time change]: Tried to change to: 2010-Oct-12, 11:00:41
79 2010-10-12, 11:03:10 [Message ]: Logbook started
80 2010-10-12, 11:06:34 [Message ]: Logbook started

```

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**ORTHOPHOS XG 3D / Ceph,
ORTHOPHOS XG 3D^{ready} / Ceph**
Compatibility list / Declaration of conformity by system integrator
 for digital radiography with SIDEXIS

System integrator (system installer)							
Company, address							
Operator							
Name, address							
Component	Requirement / Approval	X	Not X	Technical characteristics	X	Not X	Inst. regulations Function test
PC unit	Manuf.		<input type="radio"/>	IBM compatible PC with DualCore 2 GHz, RAM ≥ 2 GB, Windows XP Professional / Professional / Ultimate, ≥ 500GB free storage space		<input type="radio"/>	For the Federal Republic of Germany: On site installation according to DIN VDE 0100-710
	Type			Graphics card ≥ 1280 x 1024, ≥ 16,7 Mio colors (True Color), ≥ 512 MB Refresh rate > 70Hz		<input type="radio"/>	
	No			Operating system: according to the op. systems approved in the installation instructions.		<input type="radio"/>	
				Removable disk drive ¹ : DVD-ROM Dual Layer		<input type="radio"/>	
Monitor	Manuf.		<input type="radio"/>	Min. screen diagonal 15" with full screen, 17" with CRT		<input type="radio"/>	Display quality tested with SMPTE test image according to the enclosed instructions
	Type			Min. resolution 1280 x 1024 Refresh rate > 70Hz		<input type="radio"/>	
	No			CRT: Dot mask min. 0,28mm Full screen monitor: Pixel pitch 0,30 mm x 0,30 mm		<input type="radio"/>	
						<input type="radio"/>	
ORTHOPHOS XG 3D / Ceph, ORTHOPHOS XG 3D^{ready} / Ceph	Manuf. Sirona Dental Systems		<input checked="" type="radio"/>	see installation instructions		<input checked="" type="radio"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful
	Type					<input checked="" type="radio"/>	
	No					<input checked="" type="radio"/>	
	Manuf. Sirona Dental Systems		<input checked="" type="radio"/>	see installation instructions		<input checked="" type="radio"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful
	Type					<input checked="" type="radio"/>	
	No					<input checked="" type="radio"/>	
	Manuf. Sirona Dental Systems		<input checked="" type="radio"/>	see installation instructions		<input checked="" type="radio"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful
	Type					<input checked="" type="radio"/>	
Other devices optionally connected to the PC (Printer, LAN...)	Manuf.		<input type="radio"/>			<input type="radio"/>	Device installed according to manufacturer's instructions
	Type					<input type="radio"/>	
	No					<input type="radio"/>	
	Manuf.		<input type="radio"/>			<input type="radio"/>	Device installed according to manufacturer's instructions
	Type					<input type="radio"/>	
	No					<input type="radio"/>	
	Manuf.		<input type="radio"/>			<input type="radio"/>	Device installed according to manufacturer's instructions
	Type					<input type="radio"/>	

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Declaration acc. to Medical Devices Directive (MDD) 93/42/EEC and IEC 60601-1-1 item 3.201.4
 The undersigned herewith certifies that the above system consisting of CE marked medical devices has been assembled in accordance with the present manufacturer's instructions, that it has been tested for its intended use and that all relevant user information including pertinent information provided by the corresponding manufacturers has been provided for placing it on the market, and that this declaration shall be kept on hand for the competent authorities for a period of five years.

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sirona
The Dental Company

System integrator (system installer)								
Company, address								
Operator								
Name, address								
Component	Requirement / Approval	Notified X/yes	Technical characteristics	Notified X/yes	Inst. regulations Function test	Notified X/yes		
PC unit	Manuf.	<input type="radio"/>	IBM compatible PC with DualCore 2 GHz, RAM: ≥ 2 GB, Windows XP Professional / Ultimate, ≥ 500GB free storage space	<input type="radio"/>	For the Federal Republic of Germany: On site installation according to DIN VDE 0100-710	<input type="radio"/>		
	Type	<input type="radio"/>	Graphics card ≥ 1280 x 1024, ≥ 16,7 Mio colors (True Color), ≥ 5.12 MB Refresh rate > 70Hz	<input type="radio"/>		<input type="radio"/>		
	No	<input type="radio"/>	Operating system: according to the op. systems approved in the installation instructions.	<input type="radio"/>		<input type="radio"/>		
		<input type="radio"/>	Removable disk drive ¹ : DVD-ROM Dual Layer	<input type="radio"/>		<input type="radio"/>		
Monitor	Manuf.	<input type="radio"/>	Min. screen diagonal 15" with full screen, 17" with CRT	<input type="radio"/>	Display quality tested with SMPTE test image according to the enclosed instructions	<input type="radio"/>		
	Type	<input type="radio"/>	Min. resolution 1280 x 1024 Refresh rate > 70Hz	<input type="radio"/>		<input type="radio"/>		
	No	<input type="radio"/>	CRT: Dot mask min. 0,28mm Full screen monitor: Pixel pitch 0,30 mm x 0,30 mm	<input type="radio"/>		<input type="radio"/>		
		<input type="radio"/>		<input type="radio"/>		<input type="radio"/>		
ORTHOPHOS XG 3D / Ceph, ORTHOPHOS XG 3D ^{ready} / Ceph	Manuf.	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC	<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	Type			<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	No			<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	Manuf.	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC	<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	Type			<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	No			<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	Manuf.	Sirona Dental Systems	CE marking in compliance with Council Directive 93/42/EEC	<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
	Type			<input checked="" type="checkbox"/>	see installation instructions	<input checked="" type="checkbox"/>	For the FRG: Acceptance test according to DIN V 6885-151; for all countries: Acquisition generation with image storage and display successful	<input type="radio"/>
Other devices optionally connected to the PC (Printer, LAN...)	Manuf.		IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	Type			<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	No			<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	Manuf.		IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	Type			<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	No			<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	Manuf.		IEC 60 950-1 and CE marking according to Council Directive 89/336/EEC	<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>
	Type			<input type="radio"/>			Device installed according to manufacturer's instructions	<input type="radio"/>

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Declaration acc. to Medical Devices Directive (MDD) 93/42/EEC and IEC 60 601-1-1 Item 3.201.4
The undersigned herewith certifies that the above system consisting of CE marked medical devices has been assembled in accordance with the present manufacturer's instructions, that it has been tested for its intended use and that all relevant user information including pertinent information provided by the corresponding manufacturers has been provided for placing it on the market, and that this declaration shall be kept on hand for the competent authorities for a period of five years.

Place, date _____ Name of system integrator (in block letters) _____ Signature _____

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© Sirona Dental Systems GmbH 2012
Sirona Dental Systems GmbH